

# At Home Wound Treatment

TRUSTED BY  
*Physicians Nationwide*

## THE MEDSHIP GOAL:


To help provide innovative surgical dressings to patients by partnering with physicians. Our expertise allows us to give administrative assistance and audit support. We strive to create better outcomes, compliance, and convenience to patients.

## PATIENT GREETING:

The surgical dressings delivered to you have been prescribed by your physician. These were chosen to maximize your healing potential. It is important that you use each product as instructed.

## YOUR CARE PACKAGE:


You have received collagen dressings or collagen powder to be applied to your wound as your primary dressing. Also provided are secondary composite dressings to apply as a cover over your primary dressing. Be sure to follow your surgeon's instructions on the frequency of dressing changes.



## PRIMARY DRESSINGS:

### Collagen Dressings

Dressings with porous collagen membranes are designed to be permeable and breathable, with fluid control similar to skin. They are flexible, pliable, elastic, and durable for ease of application.



### Collagen Powder

Powder particles have a high surface area, allowing them to effectively penetrate and deliver the maximum amount of collagen into the wound surface.



If you have any questions, please contact MedShip:

**CALL:** 937-452-6115

**EMAIL:** [info@medship.health](mailto:info@medship.health)


[medship.health](https://www.medship.health)



## SECONDARY DRESSINGS:

### Silicone Composite Dressings

Secondary dressings that are waterproof, absorbent, have a non-traumatic border, and provide an ideal cover for wounds of many types. They assist in protecting the wound to prevent bacterial contamination or maceration.



### Island Composite Dressings

Secondary dressings that contain a non-adherent contact layer to cover the wound site. Consisting of a soft absorbent pad that collects exudate and protects the wound from trauma, they feature a non-woven backing with an adhesive border to gently hold the pad in place. With a waterproof backing, these dressings keep wounds dry and free of external contaminants.



[medship.health](https://www.medship.health)

# Surgical Dressing Application Instructions

The surgical dressings delivered to you have been prescribed by your physician. These were chosen to maximize your healing potential. It is important that you use each product as instructed.

You have received collagen dressings or collagen powder to be applied to your wound as your primary dressing. Also provided are secondary dressings to apply as a cover over your primary dressing. Be sure to follow your physician's instructions on the frequency of dressing changes.

## APPLICATION PROCESS:

- 1** Remove prior dressings entirely and discard.  
*Each dressing is designed for one-time use only. Do not reuse any dressing.*
- 2** Rinse the wound. Sterile saline solution or water is recommended, but warm water may be used in its place.
- 3** Pat the wound with a clean towel to remove excess fluids/water.
- 4** Apply the collagen dressing or powder as a primary dressing.  
*If you have received collagen dressings, you may cut them to best fit your wound. Both sides of the dressing are designed for coverage, so you may apply either side. If you have received collagen powder, use one packet per wound.*
- 5** Cover the primary (collagen) dressing with the secondary (silicone composite) dressing.  
*Peel back the edge of the folded release paper and center the contact layer (gauze pad) over the primary (collagen) dressing. Anchor the exposed adhesive edge to your skin while removing the remaining release paper. Gently smooth all adhesive edges onto your skin.*
- 6** Repeat the application process each subsequent day until the wound dressings are exhausted, or your wound has healed.

*Our priority is to ensure you are supported every step of the way.*



### INSTRUCTIONAL VIDEOS

Visit [medship.health/how-to-use](https://medship.health/how-to-use) to view instructional videos of how to apply your surgical dressings.

### CUSTOMER CARE

If you have any questions, please call MedShip's dedicated customer service department at 937-452-6115.



[medship.health](https://medship.health)



### PRECAUTIONS:

Do not use if you have any collagen or bovine sensitivities